

For more information, please contact :

EnvisionTEC Sales Offices

DE

EnvisionTEC GmbH
Brüsseler Straße 51
D-45968 Gladbeck
Germany
Phone +49 (0)2043-98750
<http://envisiontec.com>

U.S.A.

EnvisionTEC, Inc.
15162 S. Commerce
Dearborn, MI, 48120
USA
Phone +1-313-436-4300
<http://envisiontec.com>

Manufacturer

The Netherlands

Vertex-Dental B.V.
Centurionbaan 190
3769 AV Soesterberg
The Netherlands
Phone +31 88 616 04 40
info@vertex-dental.com
www.vertex-dental.com



<http://envisiontec.com>

E-GUIDE TINT

Monomer based on acrylic esters for manufacturing of
3D-printed surgical guides



Instruction For Use

E-Guide Tint

E-Guide Tint is a monomer based on acrylic esters for manufacturing of 3D-printed surgical guides. Suitable for printing all types of surgical guides. E-Guide Tint is a class I material and CE-certified.

The following instructions for use are for dental professionals who use E-Guide Tint as a surgical guide material. E-Guide Tint is intended exclusively for professional dental work. This instruction for use provides also information about safety and environmental aspects, a safety data-sheet is available on <http://envisiontec.com> and at local dealers. In case more information is needed about the processing of E-Guide Tint material contact the EnvisionTEC office. Also see information at the end of this document.

E-Guide Tint is a monomer based on acrylic esters for manufacturing of 3D-printed surgical guides. Suitable for printing all types of surgical guides. E-Guide Tint is a class I material and CE-certified.

IFU ID: IEGUIDETINT201703UK

Description and effects

E-Guide Tint can be used in combination with DLP based 3D printers which support EnvisionTEC materials.

Contra-indication

E-Guide Tint should not be used for any other purpose than dental work only. Any deviation from this instruction for use may have an adverse on the chemical and physical quality of E-Guide Tint. In case of an allergic reaction, please contact a medical physician.

Hazard & Precaution (H&P phrases)

Inhalation:

Irritating to respiratory system. High atmospheric concentrations may lead to irritation of the respiratory tract, dizziness, headache and anesthetic effects.

Skin contact:

May cause sensitization by skin contact. Irritating to skin, repeated and/or prolonged contact may cause dermatitis.

Eye contact:

High vapor concentration may cause irritation.

Ingestion:

Low oral toxicity, but ingestion may cause irritation of the gastrointestinal tract.

Protection:

Wear protection when handling E-Guide Tint. Protective glasses and gloves are advised. Information about the handling of the product can be found in the safety data-sheet, which is available on <http://envisiontec.com>.

Hazard Phrases:

Pictogram




Signal word Warning

Hazard statements

H315 Causes skin irritation
H317 May cause an allergic skin reaction


Processing / Post-curing

Make sure that you work as clean as possible, dirty reservoirs or machines can cause deformation and therefore failure of the printed objects!



Shake for 5 minutes before use

Before using E-Guide Tint, make sure to mix the product in the original packaging for 5 minutes. Color deviations and print failures may occur when shaken insufficiently.




Fill printer reservoir

Pour the liquid material in the reservoir of the 3D-printing machine.




For printer settings see IFU of 3D Printer

Follow the instructions for use of the printer. These are delivered together with the printer.




Remove printed parts from platform

When the machine has finished its program remove the building platform from the machine. Place the platform on some paper or cloth with the built jobs facing upwards. The printed jobs can now be removed from the platform using a suitable knife.




Cleaning Pieces step 1

Rinse the printed jobs for three minutes in an alcohol solution (>90%) to eliminate any excess material, using of an ultrasonic bath.



Cleaning Pieces step 2

Rinse for two minutes in a clean alcohol solution (>90%). Rinsing in a alcohol solution should not take longer than 5 minutes, as this may cause defects in the printed parts.



Post-Cure

After cleaning make sure the printed parts are dry and free of alcohol residual. Then place the printed jobs in a UV- light curing box for final polymerization.

Post-curing is an UV-light treatment to ensure that EnvisionTEC materials obtain full polymer conversion. Through this the residual monomer is reduced to a minimum and the required mechanical properties are obtained. This procedure is a necessary step to produce a biocompatible end-product. We strongly advice to make use of the Vertex-Dental™ LC-3DPrint Box.

Specific curing-time E-Guide Tint

EnvisionTEC Material	Time (min.)	Wavelength (nm.)	Total output Light (Watt)	UV lightbox output W*Sek=J (kJ)
E-Guide Tint	10	Blue UV-A 315-400 UV-Blue 400-550	UV-A 108 UV-Blue 108	129,6

The unit used at Vertex-Dental has 6x 18W/71 lamps (Dulux L Blue) and 6, 18W/78 lamps (Dulux blue UV-A). The calculated output is based on the UV light UVA lamp Blue. Please notice that the light sources and the printing machine need a routine maintenance following the manufacturer instructions.

Finishing

Remove any support structures and finish jobs if necessary, using conventional dental methods and instruments. Differences in color nuance may occur due to production in batches of the raw material and product, inadequate shaking of the original packaging before use or insufficient post-curing.

Storage conditions, expiry date and transport

Store the product in the original packaging at room-temperature in a dry and dark area, preferably not exceeding 25°C. Close the packaging after each use. The expiry date of the product is mentioned on the product label. In case of exceeding the expiry date, the product is no longer guaranteed in terms of treatment. Do not expose to UV-light and moisture.

Plastic and packaging waste

The product E-Guide Tint in its polymerized form is not environmentally harmful. Residual waste material in its liquid state should be delivered to a collection point for chemical waste material.

Cleaning instructions

EnvisionTEC 3D-printing material should be cleaned with non chemical products. If disinfecting before intended use is required, an ethanol solution can be used. EnvisionTEC E-Guide Tint can be sterilized by making use of an autoclave or using gamma-ray sterilisation. Do not use a dishwasher or thermal disinfectant.

EnvisionTEC recommends the following method of sterilization:

- Place the surgical guide in a standard sterilization pouch
- Apply sterilization by placing the pouch containing the guide in an autoclave for 15 minutes at 121°C or 3 minutes at 138°C. Make sure no mechanical forces are applied to the guide during sterilization.
- Let the guide cool down to room temperature before using the guide. Make sure no mechanical forces are applied to the guide during cooling down.
- Please make sure that the surgical guide is fully post-cured before sterilisation! this is important with regard to form stability.

For post-curing procedure please see § Processing / Post-curing.

Delivery units

The product E-Guide Tint is available in the following packaging size: 1000 gr.

Explanation of symbols on labelling



CE mark



Batch number of product



Manufacturer



Keep away from sunlight



Consult instructions for use



Use-by date



Temperature limit; recommended storage temperature